

# **USA LAB HeatPro™ Mantle Manual**

USA-XHM-2L, USA-XHM-5L, USA-XHM-10L, USA-XHM-20L







# **Safety Precautions and Explanations**

At USA Lab, safety is our number one priority. The following information provides guidelines for safety when using USA Lab equipment. Any piece of machinery can become dangerous to personnel when improperly operated or poorly maintained. ALL employees operating and maintaining USA Lab equipment should be familiar with its operation, thoroughly trained, and Instructed on the best safety practices.

Most industry accidents are preventable through safety awareness.

## **Training**

It is the responsibility of the customer to ensure that all personnel who will be expected to operate or maintain the equipment. Participate in training and instruction sessions to become trained operators. All personnel operating, inspecting, servicing, or cleaning this equipment must be properly trained in the operation and machine safety. **BEFORE** operating this equipment, read the operating instructions in this equipment manual. Become thoroughly familiar with the machinery and its controls.



### Safety

- Never leave the equipment running unattended and use this equipment only for its intended purpose.
- Ensure that all power sources are turned off when the machine is not in use. This encompasses electrical and pneumatic power.
- Read the manual for any special operational instructions for each piece of equipment. All
  USA Lab authored manuals are typically included with each device as well as posted online.
- Know how the equipment functions and understand the operating and halting processes.
- Wear the appropriate personal protective equipment for the task.
- When working on or around all equipment, avoid wearing loose clothing, jewelry, unrestrained long hair, loose ties, belts, scarves, or articles that may be caught in moving parts. Keep all extremities away from moving parts. Entanglement can cause death or severe injury.
- For new equipment, check input voltage and compare with the equipment voltage rating. DO NOT supply the incorrect power to any equipment for any reason whatsoever. Electrical specifications for your machine are printed on the machine tag. A properly grounded receptacle is required for safe operation regardless of voltage requirements.
- Keep the equipment operating zone free of obstacles that could cause a person to trip or
  fall toward an operating machine. Keep fingers, hands, or any part of the body out of the
  machine and away from moving parts when the machine is operating.
- Any machine with moving parts and/or electrical components can be potentially dangerous
  no matter how many safety features it contains. Stay alert and think clearly while operating
  or servicing the equipment. Be aware of operations and personnel in your surroundings. Be
  attentive to indicator lights, warning lights, and/or operator interface screens displayed on
  the machine and know how to respond.
- Do not operate machinery if you are fatigued, emotionally distressed, or under the influence of drugs or alcohol.
- Know where the FIRST AID SAFETY STATION is located.
- Know where the FIRE EXTINGUISHING EQUIPMENT is located.
- Never sit or stand on the machine or on anything that might cause you to fall against the machine.
- Rotating and moving parts are dangerous. Keep clear of the operating area. Never put any foreign object into the operating area.
- Use proper lifting and transporting devices for heavy equipment. Some types of equipment can be extremely heavy. An appropriate lifting device should be used.
- Use caution when moving portable equipment. In some cases, the machinery can be heavy and/or may be top heavy. Portable equipment can gain momentum during transport and must always be controlled.



# **Symbols and Warnings**

Below are examples of commonly used symbols and what they mean.

Understand them and their potential consequences.





#### **Product Features**

- 1. 300°C Maximum operating temperature.
- 2. PID heater control algorithm.
- 3. Auto-tuning / Auto-calibration.
- 4. Temperature probe fault detection.
- 5. Capacities from 2L-20L.
- 6. Brushless DC stirring motor.
- 7. Digital RPM display.
- 8. Over-temperature protection.
- 9. Internal and external probes.

#### Warnings

- 1. Verify that the unit is properly grounded.
- 2. Do not use heating without the stirrer function enabled.
- 3. You must use a magnetic stir bar inside the flask.
- 4. When shutting down the mantle, first disable the heating function. Leave the stirrer on until cooled, a minimum of 20 minutes.
- 5. All surfaces of the mantle can become hot. Use caution when touching any part.
- 6. On the first use, smoke and smells can come from the bowl. This is known as burn-in.
- 7. After burn-in, the bowl will become discolored. This is normal.
- 8. Perform burn-in in a well-ventilated area such as a fume hood or outdoors.
- 9. 5L, 10L, and 20L will make creaking and cracking noises when heating up. This is normal.
- 10. DO NOT USE with wet hands, in high humidity, or if something spills inside the bath.
- 11. If the unit does not power on, check the rear fuse. Replace if necessary.
- 12. If the unit is not going to be used for an extended period. Clean the unit and store it in a dry place away from corrosive gases.

# **Specifications**

SPECIFICATIONS				
Capacity	2L	5L	10L	20L
Overall Dimensions	416*282*203	380*380*275	420*420*315	494*494*315
	16.4" x 11.10" x 8"	15" x 15" x 10.8"	16.5" x 16.5" x 12.4"	19.44" x 19.44" x 12.4"
Power	110V-120V, 650W 5.9A	110V-120V, 1100W 10A	110V-120V, 2100W 19.1A	220V-240V, 3000W 13.6A
Motor	24V 40W	24V 40W	24V 40W	24V 40W
Gross Weight	22lbs.	33lbs.	48.5lbs.	55lbs.
Max Temp	300°C	300°C	300°C	300°C
Max RPM	2500RPM	2500RPM	2500RPM	2500RPM



## **Stand Assembly**

- 1. Install the vertical pole into the top of the mantle on the rear half. (left or right side)
- 2. Open and adjust the blue boss clamp. Attach loosely.
- 3. Attach finger clamp to the horizontal side of the blue boss clamp.
- 4. Adjust boss clamp and finger clamp to center of the flask neck and tighten.

## **Operation**

\*Verify that the unit is plugged in to the correct power supply voltage\*



#### **Control panel**



- A Heating On/Off
- **B** Stirring On/Off
- **C Settings / Hold for Auto-tune**
- D Shift / Move cursor
- E Lower / Reduce
- F Higher / Raise
- **G** Current temperature reading (PV)
- H Set Temperature (SV)
- I Current stirrer RPM

#### Using the unit:

- Press the heating button (A) and the stirring button (B).
- Hold **Set (**C) for 3 seconds.
- Adjust the **SV** (H) value using the **Lower** (E) and Higher (F) buttons as needed. Use **Shift** (D) to make large changes.
- Press **Set** (C) to adjust the stirrer **RPM** (I).
- Adjust the **RPM** (I) value using the **Lower** (E) and H**igher** (F) buttons as needed. Use **Shift** (D) to make large changes.
- Hold **Set** (C) for 3 seconds to set the changes.

Now the heating and stirring are set for operation. Repeat the steps above to adjust settings while running the unit. *Ramp up the temperature slowly over time.* 



## **Auto-tune / Auto-calibration**

Only perform this tuning at the desired running temperature, stirrer speed, fluid type, and fluid volume.

Hold **Shift** (D) for 8-10 seconds to enter the Auto-tune menu. Change the value to "0001" and then press **Set** (C). An "AT" symbol will flash on the screen while Auto-tune is running. The indicator will disappear when it has been completed. During the process the mantle will heat and cool in three cycles.

#### **Internal Probe**

The included accessories contain a small black plug. This part is used in the probe socket to use the internal probe instead of the external probe. To use it, unplug the external probe and replace it with the black plug.

#### **Probe Fault**

If the probe breaks an alarm will sound and the display will read "Er-1". If this happens the probe will need to be replaced.

## **Jumping Stir Bar**

If the stir bar is not stirring the fluid and is instead jumping around. Reduce the speed to a very low value. Turn the stirring function off. Wait for a minute or two. Then turn it back on. Gradually increase the speed over time.

#### **Burn-in Instructions**

This process may produce excessive smoking and discoloration of the glass fiber bowl. This is due to leftover oils from the manufacturing process.

- 1. Plug unit in, in a very well-ventilated area.
- 2. Plug in the black jumper into the probe port
  - 3. Turn on the mantle
  - 4. Set the temperature to 40°C
- 5. Once the temperature reaches the set point, increase the setting by 15°C.

(Repeat this process until 250°C is reached.)

6. Turn off the mantle and let it cool to room temperature.

### **Troubleshooting**

If you are experiencing any issues with the mantle. Please contact our Technical Department.

tech@usalab.com

or

Call (734) 855-4890 Ask for a Technician